

**ESCO**

WORLD CLASS. WORLDWIDE.



# Minimizing Contamination in Cell Culture

# Overview

- ▶ Introduction
- ▶ Different types of contamination
- ▶ Sources of contamination
- ▶ Laboratory sanitation
  - ▶ Hygiene
  - ▶ Laboratory methods and skills

# Contamination

- ▶ Contamination is a staple problem in cell culture
- ▶ All cell culture laboratories and workers have experienced contamination
- ▶ Affects both the use of cultures and research quality
- ▶ Most apparent consequences of contamination:
  - ▶ Time
  - ▶ Money
  - ▶ Effort
  - ▶ Erratic experimental results



# Consequences of contamination

## ► Loss of resources

- Time
- Money
- Effort



# Consequences of contamination

## ▶ Reproducibility crisis

### ▶ Mycoplasma contamination

- ▶ very difficult to detect via ocular inspection
- ▶ may lead to undetected contamination that may affect the current results


### ▶ Cross-contamination

- ▶ difficult to detect, unless PCR is done to confirm the identity of cells
- ▶ Cross-contaminated or misidentified cells lead to erroneous and irreproducible results.

# Types of Contamination

- ▶ Microbial Contamination

- ▶ Fungi
- ▶ Yeast
- ▶ Bacteria
  - ▶ Mycoplasma



- ▶ Cross-contamination

# Microbial Contamination

- ▶ Three most common microbial contaminants:
  - ▶ Fungi
  - ▶ Yeast
  - ▶ Bacteria
- ▶ Most common type of contamination
- ▶ Can be easily detected through microscopy by checking the medium:
  - ▶ Turbidity of the medium
  - ▶ Presence of particles in suspension
  - ▶ Rapid decline in pH (yellow color, indicates acidity)




# Mycoplasma

- ▶ Smallest bacteria
- ▶ Lack cell wall
  - ▶ Resistant to common antibiotics
  - ▶ Parasitic
- ▶ Cannot be detected through microscopy
- ▶ Presence of mycoplasma must be checked through PCR
- ▶ Not benign in cell culture; can alter the host's cell function, growth, metabolism, etc.
- ▶ Can grow to high densities without any visible signs of contamination



# Cross-contamination

- ▶ Cross-contamination of one cell line with another cell line
- ▶ Results from mishandling of cell lines during experiments and maintenance
- ▶ May lead to the replacement of the original cell with the contaminant, especially when the contaminant proliferates faster than the original cell line



# Sources of Contamination

- ▶ Indirect sources of contamination
  - ▶ Laboratory
  - ▶ Equipment
  - ▶ Personnel
- ▶ Direct sources of contamination
  - ▶ Reagents
  - ▶ Cell lines and primary cells

**Each possible source has a possible risk of transferring different types of contamination.**

# Laboratory

- ▶ Air-borne contaminants
  - ▶ Ventilation, doors, windows
- ▶ Water
  - ▶ Sink, water bath, incubator pan
    - ▶ Water microbes
    - ▶ Air-borne water microbes
- ▶ Dirt and dust
  - ▶ House germs

# Personnel

- ▶ Potential transfer of contaminants
  - ▶ Mouth and nose
    - ▶ Mycoplasmas, bacteria, viruses
  - ▶ Clothes
    - ▶ Bacteria, fungi, yeast
  - ▶ Hands
    - ▶ Bacteria, fungi, yeast
    - ▶ Mycoplasma
  - ▶ Shoes
    - ▶ Bacteria, fungi, yeast

# Reagents

- ▶ Culture medium is sterile when delivered
- ▶ FBS was one of the main sources of mycoplasma contamination
- ▶ Growth factors, peptides, antibodies could also be contamination sources
- ▶ Secondary contamination may be present in aliquots prepared in the laboratory.


# Cell lines

- ▶ Newly bought cell lines might be delivered contaminated
  - ▶ Mycoplasma and other microbes
  - ▶ Cross contamination
- ▶ Untested old cell lines might have contamination
  - ▶ Mycoplasma and other microbes
  - ▶ Cross contamination
- ▶ Primary cells may be contaminated from “semi-sterile” tissue preparation
  - ▶ Mycoplasma and other microbes

# Laboratory sanitation

- ▶ Indirect sanitation
  - ▶ Laboratory and equipment
  - ▶ Personnel
- ▶ Direct sanitation
  - ▶ Reagents
  - ▶ Cell lines and primary cells

# Cell culture facilities and labs must be sterile.





# Cell culture laboratory

- ▶ Dedicated air filtration system
  - ▶ Positive pressure to keep microbes out of the lab
- ▶ PPE
  - ▶ Lab gown, gloves, shoes
  - ▶ Gowning rooms
- ▶ High cleaning frequency
- ▶ Sufficient sanitizers and disinfectants

# CO<sub>2</sub> Incubators

- ▶ Fast recovery time after door openings
- ▶ Sterile surface
- ▶ Minimize door openings
- ▶ Inner doors
  - ▶ Reduce entry of microbes
  - ▶ Keep pH, temperature, humidity constant



# Cleaning and Maintenance of CO<sub>2</sub> Incubators

- ▶ Decontaminate incubators at least once a month
- ▶ Clean the incubator once a month
  - ▶ Disassemble and clean the inner parts
  - ▶ Clean the inner walls
  - ▶ Rinse parts with distilled water
- ▶ Consider the use of copper chamber
  - ▶ Inhibits bacterial growth better than stainless steel



# Humidity pan – CO2 Incubators


- ▶ Air-borne microbes might proliferate in the humidity pan
- ▶ Fan-less design can minimize the distribution of microbes in the chamber
- ▶ Autoclaved water must be used
- ▶ Antibiotics may be added to the water to inhibit germ growth



# Water bath

- ▶ Used for thawing cells, warming of media and other reagents
- ▶ A significant hub of microbes
  - ▶ Ensure that the water bath has lid
  - ▶ Clean regularly – once a week
  - ▶ Use supplements such as antibiotics to prevent growth of germs
  - ▶ Disinfect flasks, tubes, bottles well after removal from water bath and before using for aseptic techniques

# Pipettes and pipette tips



- ▶ Designate pipettes to be used in the cell culture lab only
- ▶ Wipe pipettes before using inside the biological safety cabinet
- ▶ Disassemble and clean the pipettes regularly to prevent germ growth
- ▶ Check for the accuracy of the pipettes through calibration at least once every year
- ▶ Dispose pipette tips properly
  - ▶ Do not leave used tips inside the biological safety cabinet

# Other lab equipment

- ▶ Wipe off surfaces daily if possible
- ▶ Check if the lids of equipment not in used are closed
- ▶ Immediately clean all spills
  - ▶ Autoclave rotors of centrifuges



**ESCO**

WORLD CLASS. WORLDWIDE.

# Personnel


- ▶ Wash hands thoroughly before entry into the cell culture lab
- ▶ Ensure PPE (Personal Protection Equipment) immediately
  - ▶ Lab gown
  - ▶ Gloves
  - ▶ Shoes
  - ▶ Mask (if needed, depends on BSC design)
- ▶ Wash hands before leaving the cell culture lab





# Cell lines and primary cells

- ▶ Acquire cell lines only from authenticated and reliable sources
  - ▶ European Collection of Authenticated Cell Cultures (ECACC)
  - ▶ American Type Culture Collection (ATCC)
- ▶ Quarantine cells when testing
  - ▶ Bacteria, fungi, yeast (in culture)
  - ▶ Mycoplasma (PCR)
  - ▶ Cross-contamination (STR service for human cell lines)
- ▶ Primary cells must be maintained in a separate incubator



# Biological safety cabinets

- ▶ Regular cleaning and maintenance
- ▶ Disinfect the surface before use
- ▶ Disinfect materials before transfer
- ▶ Arrange materials inside the BSC strategically
  - ▶ Only needed materials should be placed inside the BSC
  - ▶ Enough working space must be maintained
- ▶ Spray disinfectant after use



# Proper pipetting

## Aspirate liquid

- ▶ Vertical hold
- ▶ Aspirate slowly for precision and to avoid pipette contamination

## Release liquid

- ▶ Maintain a 45°C for sterility
- ▶ Release at medium speed to avoid droplets

# Filtered micropipette tips

- ▶ Filter tips for increased protection
  - ▶ Protects both pipette and culture from contamination
  - ▶ Use particularly for important cell culture steps, such as cell banking



# Vessels



- ▶ Lids of flasks can be an easy contamination route
  - ▶ Place lids safely when working outside the BSC
  - ▶ Upside up or upside down both allowed inside the BSC when done properly

# Labeling

- ▶ Label flasks, dishes, plates, and tubes clearly
  - ▶ Prevents cross-contamination
- ▶ Ethanol-resistant pen must be used
  - ▶ Prevents loss of label when disinfecting

# Safe transport

- ▶ Slightly tilt the culture flasks upwards to avoid medium contact
- ▶ Avoid spillage when transporting plates and dishes

Prevent contamination before  
it's too late!

**ESCO**

WORLD CLASS. WORLDWIDE.





**ESCO**  
WORLD CLASS. WORLDWIDE.  
[www.escoglobal.com](http://www.escoglobal.com)

Life science tools and equipment  
for your cell culture laboratory.